

ExCell Bio

OptiVitro[®] NK Cell Expansion Serum-free Basic Kit P01

For Research and Manufacturing Use
Not Intended for Diagnostic and Therapeutic Use

User Manual

Catalog Number NE000-N032
NE000-N031
NE000-N031S



Product description

OptiVibro® NK Cell Expansion Serum-free Basic Kit P01 has been specifically designed for the in vitro expansion of human Natural Killer (NK) cells derived from either peripheral blood mononuclear cells (PBMCs) or cord blood mononuclear cells (CB-MNCs). The kit is composed of three main components: OptiVibro® NK Cell Basal SF Medium P01, OptiVibro® NK Cell SF Medium Supplement, and OptiVibro® Cytokine III. All of these components are serum-free, xeno-free, and have been manufactured in strict compliance with GMP regulations. This Kit can be combined with OptiVibro® NK Cell Expansion Serum-free Kit P01 (NE000-N02#) to facilitate and support the expansion of NK cells in vitro.

Contents and storage

Catalog No.	Product name	Amount	Storage	Shelf life ^[1]
NE000-N032	OptiVibro® NK Cell Expansion Serum-free Basic Kit P01	1 kit	-	-
BA0092	OptiVibro® NK Cell Basal SF Medium P01	1000 mL	2-8 °C Protect from light	12 months
BA0102	OptiVibro® NK Cell SF Medium Supplement	8 mL	2-8 °C Protect from light	18 months
BA0132	OptiVibro® Cytokine III	310 µL	-20 °C	12 months
NE000-N031	OptiVibro® NK Cell Expansion Serum-free Basic Kit P01	1 kit	-	-
BA0091	OptiVibro® NK Cell Basal SF Medium P01	500 mL	2-8 °C Protect from light	12 months
BA0101	OptiVibro® NK Cell SF Medium Supplement	4 mL	2-8 °C Protect from light	18 months
BA0131	OptiVibro® Cytokine III	155 µL	-20 °C	12 months
NE000-N031S	OptiVibro® NK Cell Expansion Kit Basic P01 (Sample)	1 kit	-	-
BA0091S	OptiVibro® NK Cell Basal SF Medium P01 (Sample)	100 mL	2-8 °C Protect from light	12 months
BA0101S	OptiVibro® NK Cell SF Medium Supplement (Sample)	0.8 mL	2-8 °C Protect from light	18 months
BA0131S	OptiVibro® Cytokine III (Sample)	31 µL	-20 °C	12 months

^[1] The Shelf-Life may be extended after strict validation by QC.

| Instructions for use

Prepare media

1. Place OptiVibro[®] NK Cell Basal SF Medium P01 and OptiVibro[®] NK Cell SF Medium Supplement under a sterile laminar flow hood.
2. Add 4 mL/8 mL OptiVibro[®] NK Cell SF Medium Supplement to 500 mL/1000 mL OptiVibro[®] NK Cell Basal SF Medium P01.
3. Tighten and mix the complete OptiVibro[®] NK Cell Basal SF Medium P01 thoroughly.

Note: It is recommended to use complete OptiVibro[®] NK Cell Basal SF Medium P01 within four weeks after mixed.

4. To prepare the complete medium (short name: "NK-SFM"), add 155 μ L/310 μ L of OptiVibro[®] Cytokine III to 500 mL/1000 mL of the previously prepared OptiVibro[®] NK Cell Expansion Serum-free Medium P01.

Note:

1) The complete medium (the short name is 'NK-SFM' in the following protocol) is stable for 3 weeks when stored at 2-8°C in the dark.

2) OptiVibro[®] Cytokine III can be aliquoted for small volume culture use but should be limited to three freeze-thaw cycles.

Culture NK cells from PBMCs

OptiVibro[®] NK Cell Expansion Serum-free Kit P01 is designed for culturing NK cells from peripheral blood mononuclear cells (PBMCs), cord blood mononuclear cells (CB-MNCs), or NK cells derived from iPS cells. It can be used together with OptiVibro[®] NK Cell Expansion Serum-free Kit P01 (NE000-N02#) for feeder-free NK cell culture systems. The kit does not necessarily require serum or serum replacement, but supplementing with heat-inactivated autologous plasma, serum replacement, or human AB serum can increase cell expansion folds. This protocol outlines the procedures for culturing NK cells from PBMCs, starting with a T75 flask as an example.

1. Thaw OptiVibro[®] Cytokine I (OptiVibro[®] NK Cell Expansion Serum-free Kit P01, NE000-N02#) at room temperature one day before NK cell activation. Add 45 μ L OptiVibro[®] Cytokine I to 15 mL sterile DPBS and mix well. Transfer the mixed liquid to a sterile T75 flask, shaking it slightly to ensure the liquid covers the bottom of the flask. Store the flask at 2-8°C overnight.
2. Prepare fresh PBMCs following standard PBMC separation protocols or quickly thaw (<1 minute) frozen vials of PBMC cells in a 37°C water bath.

3. If using fresh PBMCs, wash them with sterile DPBS and use them directly. If using frozen cells, thaw them one day before NK cell activation, place them at a concentration of around 2×10^6 cells/mL in complete OptiVibro® NK Cell Expansion Serum-free Medium P01 without extra cytokines, and incubate them in a humidified 37°C incubator with an atmosphere of 5% CO_2 for 16-24 h.
4. It is optional to sort NK cells using magnetic beads with antibodies before the activation.
5. Centrifuge cells at $400 \times g$ for 10 minutes and discard the supernatant.
6. Equilibrate the T75 flask coated with OptiVibro® Cytokine I (prepared in step 1) at room temperature and remove the liquid.
7. Equilibrate complete OptiVibro® NK Cell Expansion Serum-free Medium P01 to room temperature before use. Resuspend the PBMCs at a concentration of $2.0\text{-}2.5 \times 10^6$ cells/mL in 15 mL of complete OptiVibro® NK Cell Expansion Serum-free Medium P01 with OptiVibro® Cytokine III (NK-SFM) supplemented with 10% heat-inactivated autologous plasma.
8. Transfer the cells (from step 7) to the T75 flask (from step 6), add 150 μL of OptiVibro® Cytokine II to the medium, and shake it slightly. Incubate the cells in a humidified 37°C incubator with an atmosphere of 5% CO_2 .
9. On Day 3 after NK cell activation, feed the cells with 15 mL of NK-SFM supplemented with 10% heat-inactivated autologous plasma.
10. On Day 5 after NK cell activation, feed the cells and adjust the cell concentration to $1.0\text{-}1.5 \times 10^6$ cells/mL with NK-SFM supplemented with 5% heat-inactivated autologous plasma.
11. From Day 7 after NK cell activation, feed the cells and adjust the cell concentration to $1.0\text{-}1.5 \times 10^6$ cells/mL with NK-SFM supplemented with 1% heat-inactivated autologous plasma every 2-3 days. The cells can be transferred to bioreactors for further expansion at around Day 9-11 after NK cell activation.

Note: If culturing NK cells from CB-MNCs, it is recommended to use freshly prepared CB-MNCs instead of frozen ones to achieve higher expansion folds.